elite® P&P Regular Body

Normal Set

A-silicones Precision & Performance





USE

- Fixed prostheses
- Recommended for the Single-stage impression technique or monophase with individual tray

CHARACTERISTICS

- 1:1 Mixing ratio
- Hydrophilic
- Thixotropic
- Ideal working time/setting time ratio (snap-set effect):
 - Long working time;
 - Shorter permanence in the oral cavity
- Biocompatible
- Purple colour

ADVANTAGES

- Maximum precision in detail reproduction
- · High tear strength
- Excellent elastic recovery
- High dimensional stability
- Safe disinfection with Zeta 7 Solution and Zeta 7 Spray
- 3-year shelf life

Technical data	Normal Set
ISO 4823	Type 2 Medium-Bodied
ADA 19	Consistency
Mixing ratio Base:Catalyst	1:1
Mixing time	Tube Version: 0:30 min/ automix
Total working time (23°C/73°F)	2:00 min
Time in mouth	3:30 min
Setting time (23°C/73°F)	5:30 min
Strain in compression (min - max)	3 - 5 %
Elastic recovery	> 99.5 %
Dimensional stability (after 24 h)	< 0.20 %
Hardness (Shore A)	45 ± 2

Equipment	
C303230	Alghamix II - mechanical mixer for alginates, plasters and silicones 230V
C303215	Alghamix II - mechanical mixer for alginates, plasters and silicones 115V
C303231	Alghamix II - mechanical mixer for alginates, plasters and silicones 230V + UK plug

Packaging	
C206021	

C206021	elite® P&P Regular Normal Set 2 x 60 ml tubes (Base + Catalyst)
C206012	elite® P&P Regular Normal Set 2 x 50 ml cartridges (Base + Catalyst) + 6 mixing tips small/yellow

Accessories

C207202	Mixing pad (12 sheets)
C202100	D2 Dispenser -1:1
C207200	Spatula for silicones

KIT	
C206100	elite® P&P Intro Kit Light Normal Set - 2 x 250 ml tubs Putty Soft Normal (Base + Catalyst) + 2 x 60 ml tubes Light Body Normal + Mixing pad
C206101	elite® P&P Intro Kit Regular Normal Set - 2 x 250 ml tubs Putty Soft Normal (Base + Catalyst) + 2 x 60 ml tubes Regular Body Normal + Mixing pad
C206102	elite® P&P Intro Kit Light Normal Set - 2 x 250 ml tubs Putty Soft Normal (Base + Catalyst) + 2 x 50 ml cartridges Light Body Normal + D2 Dispenser + 6 mixing tips small/yellow - 2 Oral mixing tips
C206111	elite® P&P Trial Kit Regular Normal Set - 2 x 100 ml tubs Putty Soft Normal + 2 x 60 ml tubes Regular Body Normal