

elite® P&P Light Body

Normal & Fast Set

A-silicones Precision & Performance



USE

- Fixed prostheses
- Suggested techniques: two-stage or single-stage technique with two viscosities

CHARACTERISTICS

- 1:1 Mixing ratio
- Hydrophilic
- Thixotropic
- Ideal working time/setting time ratio (snap-set effect):
 - Long working time;
 - Shorter permanence in the oral cavity
- Biocompatible
- Colour:
 - Light Blue (Normal Setting)
 - Orange (Fast Setting)

ADVANTAGES

- Maximum precision in detail reproduction
- High tear strength
- Excellent elastic recovery
- High dimensional stability
- Safe disinfection with Zeta 7 Solution and Zeta 7 Spray
- 3-year shelf life

Technical data

	Normal Set	Fast Set
ISO 4823	Type 3 Light-Bodied	Type 3 Light-Bodied
ADA 19	Consistency	Consistency
Mixing ratio Base:Catalyst	1:1	1:1
Mixing time	Tube version: 0:30 min/ automix	automix
Total working time (23°C/73°F)	2:00 min	1:30 min
Time in mouth	3:30 min	2:30 min
Setting time (23°C/73°F)	5:30 min	4:00 min
Strain in compression (min - max)	3 - 5 %	3 - 5 %
Elastic recovery	> 99.5 %	> 99.5 %
Dimensional stability (after 24 h)	< 0.20 %	< 0.20 %
Hardness (Shore A)	45 ± 2	45 ± 2

Equipment

C303230	Alghamix II - mechanical mixer for alginates, plasters and silicones 230V
C303215	Alghamix II - mechanical mixer for alginates, plasters and silicones 115V
C303231	Alghamix II - mechanical mixer for alginates, plasters and silicones 230V + UK plug

Packaging

C206020	elite® P&P Light Normal Set 2 x 60 ml tubes (Base + Catalyst)
C206010	elite® P&P Light Normal Set 2 x 50 ml cartridges (Base + Catalyst) + 6 mixing tips small/yellow
C206011	elite® P&P Light Fast Set 2 x 50 ml cartridges (Base + Catalyst) + 6 mixing tips small/yellow

Kit

C206100	elite® P&P Intro Kit Light Normal Set - 2 x 250 ml tubs Putty Soft Normal (Base + Catalyst) + 2 x 60 ml tubes Light Body Normal + Mixing pad
C206101	elite® P&P Intro Kit Regular Normal Set - 2 x 250 ml tubs Putty Soft Normal (Base + Catalyst) + 2 x 60 ml tubes Regular Body Normal + Mixing pad
C206102	elite® P&P Intro Kit Light Normal Set - 2 x 250 ml tubs Putty Soft Normal (Base + Catalyst) + 2 x 50 ml cartridges Light Body Normal + D2 Dispenser + 6 mixing tips small/yellow - 2 Oral mixing tips
C206110	elite® P&P Trial Kit Light Normal Set - 2 x 100 ml tubs Putty Soft Normal + 2 x 60 ml tubes Light Body Normal

Accessories

C207202	Mixing pad (12 sheets)
C202100	D2 Dispenser -1:1
C207200	Spatula for silicones
C207110	Syringe Zhermack for elastomers